


# Dossier d'information Euro Pharmat

## DISPOSITIF MEDICAL


### SERINGUES TROIS PIECES LUER LOCK

<b>1. Renseignements administratifs concernant l'entreprise</b>		<b>Date de mise à jour : Mars 2016</b>
1.1	Nom : <b>TERUMO France</b>	
1.2	Adresse complète : Bât.Renaissance 3 Rond-Point des Saules 78284 Guyancourt Cedex	Tel: <b>01 30 96 13 00</b> Fax : <b>01 30 43 60 85</b> e-mail : <b>terumo.france@terumo-europe.com</b> Site internet : <b>www.terumo-europe.com</b>
1.3	Coordonnées du correspondant matériovigilance : Sara DELANNAY	Tel : <b>01 30 96 13 03</b> Fax : <b>01 30 43 60 85</b> e-mail : <b>sara.delannay@terumo-europe.com</b>
<b>2. Informations sur le dispositif ou équipement</b>		
2.1	<b>Dénomination commune</b> : selon la nomenclature d'Europharmat® <b>Seringue</b>	
2.2	<b>Dénomination commerciale</b> : <b>Seringue trois pièces Luer Lock</b>	
2.3	<b>Code nomenclature</b> : <b>Code GMDN</b> : seringue : 35904 et seringue opaque : 45492 <b>Code CLADIMED</b> : <b>K54BB03</b>	
2.4	<b>Code LPP*</b> : <b>Non applicable</b> * « liste des produits et prestations remboursables » inscrits sur la liste prévue à l'article L 165-1	
2.5	<b>Classe du DM</b> : <b>IIa</b> Directive de l'UE applicable : <b>93/42/CE</b> Selon Annexe n° <b>I et II</b> Numéro de l'organisme notifié : <b>CE 0197 (TÜV Rheinland, Cologne Allemagne)</b> Date de première mise sur le marché dans l'UE : <b>Avant 1998</b> Fabricant du DM : <b>Terumo Belgique, Terumo Japon et Terumo Philippines</b>	
2.6	<b>Descriptif du dispositif (avec photo, schéma, dimensions, volume, ...)</b> :  Seringues trois pièces, stériles et non pyrogènes, conformes à la norme NF EN ISO 7886-1 et -2 (seringue/pousse seringue)  Elles sont constituées : <ul style="list-style-type: none"> <li>○ d'un corps transparent doté d'un embout luer lock centré et d'un bourrelet d'arrêt interne en haut du corps pour éviter la sortie du piston en bout de course. Graduation conforme à la norme, impression noire ou bleue.</li> <li>○ d'un piston muni d'un joint à double lèvre, assurant une étanchéité parfaite l'épaisseur optimale du joint, permet une mobilité du piston sans effort et sans à coup.</li> <li>○ L'intérieur du corps de la seringue et le joint sont siliconés.</li> </ul>	
		

### 2.7 Références Catalogue :

**Tableau des références**

Références	Description	Nbre Unit/boîte	Nbre unit/carton
1SS02LE1	2.5ml- Luer Lock	100	2400
1SS05LE1	5.0ml- Luer Lock	100	1600
1SS10LE1	10ml- Luer Lock	100	800
1SS30LE1	30ml- Luer Lock	50	400
8SS03L1	3.0ml -Luer lock	100	1800
8SS05L1	5.0ml- Luer lock	100	1200
8SS10L1	10ml- Luer lock	100	1200
8SS20L1	20ml- Luer lock	50	600
8SS30L1	30ml- Luer lock	25	200
8SS50L1	50ml- Luer lock	25	100
8SS50LB1	50ml - opaque	25	100



#### Conditionnement/Emballages :

UCD (Unité de commande): **100, 50 ou 25 seringues selon le volume**

CDT (Multiple de l'UCD): **Quantité variable selon le volume voir tableau ci dessus**

QML (Quantité minimale de livraison) : **Le carton**

Code à barres : **EAN 128**

#### Descriptif de la référence

POSITIONS	REFERENCES	EXPLICATIONS
1	1, 2, 8	Lieu de fabrication : 1=Japon, 2=Belgique et 8=Philippines
2-3	SS/BS	Seringue
4-5	02, 03, 05, 10, 20, 30, 50	Volume en ml : 02=2,5 03=3ml 05=5ml 10=10ml 20=20ml 30=30ml 50=50ml
6	L	Luer Lock
7-8	E1 ou 1	Stérilisation faisceau électrons

**Etiquetage** : Voir ANNEXES

**2.8** Composition du dispositif et Accessoires : pour chaque élément ou composant, précisé :

<u>Dispositif</u>	<u>Eléments</u>	<u>Matériaux</u>
<b>Seringue</b>	Corps	Polypropylène
	Piston	Polypropylène
	Joint	Elastomère thermoplastique,
	Lubrifiant	Huile de silicone

Silicone : Conforme à la pharmacopée européenne

Pour les composants susceptibles d'entrer en contact avec le patient et/ou les produits administrés, précisions complémentaires :

- Absence de Latex (Annexe 7)
- Absence de produit d'origine animale ou biologique
- Absence de PVC/Phtalates

Toutes mentions jugées utiles pour les précautions d'utilisation :

- Vérifier l'intégrité du protecteur individuel de stérilité avant utilisation
- Strict usage unique, détruire après usage selon les procédures locales d'élimination des déchets de soins

Domaine - Indications :

**2.9** Domaine d'utilisation (selon liste Europharmat) : Injection

Indications (selon liste Europharmat) : Injection manuelle ou pousse seringue selon références



# Dossier d'information Euro Pharmat

## DISPOSITIF MEDICAL


<b>3. Procédé de stérilisation :</b>	
	<b>DM stérile : OUI</b> <b>Mode de stérilisation du dispositif :</b> Références avec un 2 devant = Stérilisé à l'Oxyde d'éthylène, validation du process selon la norme EN ISO 11135-1:2007 Références avec un 1 ou un 8 devant= Par faisceau d'électrons Selon la norme EN ISO 11137-1:2006/11137-2:2007
<b>4. Conditions de conservation et de stockage</b>	
	Conditions normales de conservation & de stockage précautions particulières: - <b>Eviter le stockage à des températures excessives et à l'humidité</b> Durée de la validité du produit: <b>5 ans</b> Présence d'indicateurs de température s'il y a lieu: <b>Non</b>
<b>5. Sécurité d'utilisation</b>	
<b>5.1</b>	<b>Sécurité technique : Voir Annexe 1</b>
<b>5.2</b>	<b>Sécurité biologique (s'il y a lieu) : Non applicable</b>
<b>6. Conseils d'utilisation</b>	
<b>6.1</b>	<b>Mode d'emploi : Voir annexe 1</b>
<b>6.2</b>	<b>Indications : Voir annexe 1</b>
<b>6.3</b>	<b>Précautions d'emploi : Voir Annexe 1</b>
<b>6.4</b>	<b>Contre- Indications : Voir Annexe 1</b>
<b>7. Informations complémentaires sur le produit</b>	
	<b><u>Bibliographie, rapport d'essais cliniques, ou d'études pharmaco-économiques, amélioration du service rendu : recommandations particulières d'utilisation (restrictions de prise en charge, plateau technique, qualification de l'opérateur, etc) ... :</u></b>
<b>8. Liste des annexes au dossier (s'il y a lieu)</b>	
	<ul style="list-style-type: none"><li>✓ Précautions d'emploi sur la boîte (Annexe 1)</li><li>✓ Boîte, Etiquetage blister (Annexes 2, 3)</li><li>✓ Certificat de marquage CE Japon (Annexe 4)</li><li>✓ Certificat de marquage CE Philippines (Annexe 5)</li><li>✓ Déclaration de conformité Japon (Annexe 7)</li><li>✓ Déclaration de conformité Philippines (Annexe 6)</li><li>✓</li></ul>
<b>9. Images (s'il y a lieu)</b>	
	Format gif, jpeg, png

### ANNEXE 1


#### Précautions d'emploi

**STERILE**

Sterilized by electron beam / Sterilisé par faisceau d'électrons / Sterilisiert durch Bestrahlung / Esterilizado por haz de electrones / Sterilizzato a fascio di elettroni / Esteriliseerd door electron beam / Elektronstrålesterilisert



For single use only / Strict usage unique / Nur zum sinnigen Gebrauch / Valido para un solo uso / Monouso / Voor éénmalig gebruik / Endast för enångsbruk



LATEX FREE / SANS LATEX / LATEXFREI / SIN LATEX / PRIVO DI LATTICE / LATEX-VRIJ / LATEX FRI

•Sterile and non-pyrogenic in an unopened and undamaged unit package.  
**<PRECAUTIONS>** •Do not use if the unit package or the product has been damaged or soiled.  
 •Use immediately after opening the unit packaging. •Dispose of safely after single use to avoid risk of infection. •Do not use for high pressure injection of contrast media. •**DO NOT STORE AT EXTREME TEMPERATURES AND HUMIDITY, AVOID DIRECT SUNLIGHT.** **<INSTRUCTIONS FOR USE>** To open the blister package, peel back the top layer starting from the arrow. (Fig.1)

•Steril et apyrogène dans un emballage individuel non ouvert et non endommagé.  
**<PRECAUTIONS>** •Ne pas utiliser si l'emballage individuel ou le produit a été endommagé ou souillé. •Utiliser immédiatement après ouverture de l'emballage individuel. •Éliminer de façon appropriée après usage unique pour éviter le risque d'infection. •Ne pas utiliser pour l'injection de liquide de contraste sous forte pression. •**ÉVITER LE STOCKAGE A DES TEMPERATURES EXTREMES ET A L'HUMIDITÉ. ÉVITER LA LUMIÈRE DIRECTE DU SOLEIL.** **<MODE D'EMPLOI>** Pour ouvrir le blister, détacher les deux parties, en partant de la flèche. (Fig.1)

•Steril und pyrogenfrei in ungeöffneter und unbeschädigter Einzelverpackung.  
**<VORSICHTSMAßNAHMEN>** •Nicht verwenden, wenn die Einzelverpackung oder das Produkt beschädigt oder verschmutzt sind. •Nach Öffnen der Einzelverpackung das Produkt umgehend verwenden. •Nach einmaligem Gebrauch sicher entsorgen um Infektionsrisiko zu vermeiden. •Nicht geeignet für Druckinjektionen von Kontrastmitteln. •**VERMEIDEN SIE EXTREME TEMPERATUREN UND FEUCHTIGKEIT WÄHREND DER LAGERUNG. VOR DIREKTER SONNENBESTRAHLUNG SCHÜTZEN.** **<GEBRAUCHSANLEITUNG>** Öffnen der Blisterverpackung: Ziehen Sie die obere Schicht beim Pfeil beginnend ab. (Fig.1)

•Estéril y apirógeno si el envase unitario no ha sido abierto ni deteriorado.  
**<PRECAUCIONES>** •No utilizar si el envase unitario o el producto están manchados o dañados. •Utilizar inmediatamente después de abrir el envase unitario. •Usar una vez y destruir. El uso compartido conlleva riesgo de infección. •No utilizar para inyectar medio de contraste a altas presiones. •**NO ALMACENAR A TEMPERATURAS EXTREMAS NI EN LUGARES HÚMEDOS. EVITAR LA LUZ SOLAR DIRECTA.** **<INSTRUCCIONES DE USO>** Para abrir el envase blister: tirar de las lengüetas hacia fuera siguiendo la indicación de la flecha. (Fig.1)

•Sterile og apyrogenic i en uåbnet og uskadedt styckförpackning.  
**<FÖRSIKTIGHETSÅTGÄRDER>** •Får ej användas om styckförpackning eller produkt har skadats eller blivit smutsad. •Produkten används omedelbart efter det att styckförpackningen har öppnats. •Efter engångsanvändning, kassera på ett säkert och ansvarsfullt sätt för att undvika infektionsrisk. •Får ej användas för högttrycksinjektion av kontrastmedla. •**FÖRVARAS EJ VID EXTREM TEMPERATUR ELLER LUFTHUIGTIGHET. UNDVIK DIREKT SOLLJUS.** **<BRUKSANVISNING>** För att öppna blisterförpackningen: drag iår den övre delen av förpackningen med start från pilen. (Fig.1)




Fig.1

**TERUMO CORPORATION** TOKYO 151-0072, JAPAN MADE IN JAPAN

**EC REP** TERUMO EUROPE N.V., Interleuvenaan 40, 3001 LEUVEN, BELGIUM



# Dossier d'information Euro Pharmat

## DISPOSITIF MEDICAL

### ANNEXE 2

#### Boite



# Dossier d'information Euro Pharmat


## DISPOSITIF MEDICAL

### ANNEXE 3

#### Etiquetage blister



*N° Lot et date de péremption sur le film du blister*



TÜVRheinland®

**EC Certificate**  
Directive 93/42/EEC Annex II, excluding Section 4  
Full Quality Assurance System  
Medical Devices

**Registration No.:** HD 60077473 0001

**Report No.:** 12018187 001

**Manufacturer:** Terumo Corporation  
44-1, 2-chome, Hatagaya  
SHIBUYA-KU, TOKYO 151-0072  
JAPAN


**Products:** see attachement for products included  
Replaces Approval, Registration No.: HD 60026344 0001


**Expiry Date:** 2017-08-29

The Notified Body hereby declares that the requirements of Annex II, excluding section 4 of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex II, section 5 of the aforementioned directive. For placing on the market of class III devices covered by this certificate an EC design-examination certificate according to Annex II, section 4 is required.

**Effective Date:** 2013-05-31

**Date:** 2013-05-31



Notified Body  
  
Dr. H. Lüdemann

**TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg**

TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

10.020.04.00 © TÜV, TÜV and TÜV are registered trademarks. Use without approval requires prior approval.





Doc. 1/2, Rev .0

**TÜV Rheinland**  
**LGA Products GmbH**  
**Tillystraße 2, 90431 Nürnberg**

**Attachment to  
Certificate**

**Registration No.:** HD 60077473 0001  
**Report No.:** 12018187 001

**Manufacturer:** Terumo Corporation  
44-1, 2-chome, Hatagaya  
SHIBUYA-KU, TOKYO 151-0072  
JAPAN

**Products included:**

- Blood Bags
- Blood Donor Set with/without Blood Transfusion Filter
- Blood Transfusion Filter
  
- Intravenous Catheter
- Intravenous Administration Set
- Hypodermic Syringe
- Winged Needle
- Dental Needle
- Spinal Needle
- Other Medical Needle
- Blood Administration Set
- Lancet

**Date:** 2013-05-31



**Notified Body**

  
**Dr. H. Lüdemann**



Doc. 2/2, Rev .0

**TÜV Rheinland**  
**LGA Products GmbH**  
**Tillystraße 2, 90431 Nürnberg**

**Attachment to  
Certificate**

**Registration No.:** HD 60077473 0001  
**Report No.:** 12018187 001

**Manufacturer:** Terumo Corporation  
44-1, 2-chome, Hatagaya  
SHIBUYA-KU, TOKYO 151-0072  
JAPAN

**Products included:**

- Extra-corporeal Membrane Oxygenator
- Cardiopulmonary Bypass Arterial Line Blood Filter
- Heart-Lung Bypass Defoamer
- Cardiotomy Reservoir
- Cardiopulmonary Bypass Blood Reservoir
- Haemoconcentration Filter
- Centrifugal Pump
- Angiographic Catheter
- Balloon Dilatation Catheter
- Catheter Guide Wire
- Guiding Catheter
- Catheter Introducer
- Stents
- Extension Tube
- Temperature Control Unit for Heart-Lung Bypass System  
Module
- Infusion Pump
- Syringe Infusion Pump
- Clinical Electronic Thermometer
- Clinical Electronic Blood-Pressure Monitor
- Endoscopic Electromechanical Surgical Systems

**Date:** 2013-05-31



**Notified Body**

  
**Dr. H. Lüdemann**



**EC Certificate**  
**Directive 93/42/EEC Annex V**  
**Production Quality Assurance**  
**Medical Devices**

**Registration No.:** DD 60108472 0001

**Report No.:** 12031276 001

**Manufacturer:** Terumo (Philippines) Corporation  
124 East Main Avenue  
Laguna Technopark, Binan,  
Laguna, 4024  
Philippines

**Products:** See attachments for products and sites included  
Replaces Approval, Registration No.: DD 60083914 0001

**Expiry Date:** 2021-02-11

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

**Effective Date:** 2016-02-12

**Date:** 2016-02-12



**TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg**

TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



Doc. 1/2, Rev. 0

**TÜV Rheinland**  
**LGA Products GmbH**  
**Tillystraße 2, 90431 Nürnberg**

**Attachment to  
Certificate**

**Registration No.:** DD 60108472 0001  
**Report No.:** 12031276 001

**Manufacturer:** Terumo (Philippines) Corporation  
124 East Main Avenue  
Laguna Technopark, Binan,  
Laguna, 4024  
Philippines


**Products included:**

- Syringes with Needles
- Intravenous Catheters
- Safety Needles
- Syringes with Safety Needles
- Syringes without Needles
- Hypodermic Needles

**Aspects of manufacturing concerned with securing and  
maintaining sterile conditions:**

- Urinary Drainage Bags
- Syringes for Oral / Enteral

**Date: 2016-02-12**



**Dipl.-Ing. S. Pane**



Doc. 2/2, Rev. 0

**TÜV Rheinland**  
**LGA Products GmbH**  
**Tillystraße 2, 90431 Nürnberg**

**Attachment to  
Certificate**

**Registration No.:** DD 60108472 0001  
**Report No.:** 12031276 001

**Manufacturer:** Terumo (Philippines) Corporation  
124 East Main Avenue  
Laguna Technopark, Binan,  
Laguna, 4024  
Philippines

**Manufacturing site included:**

Terumo (Philippines) Corporation  
128 East Main Avenue, Laguna Technopark, Binan, Laguna,  
4024, Philippines  
- Intravenous Catheter  
- Safety Needles  
- Syringes with Safety Needles

Aspects of manufacturing concerned with securing and  
maintaining sterile conditions:  
- Urinary Drainage Bags

**Sterilization (Electron Beam Irradiation) site included:**

Terumo (Philippines) Corporation  
124 East Main Avenue, Laguna Technopark, Binan, Laguna,  
4024, Philippines

**Date:** 2016-02-12





### **TERUMO (PHILIPPINES) CORPORATION**

124 East Main Ave., Laguna Technopark, Biñan, Laguna, Philippines  
Tel. No. (049) 541-2111 • Fax No. (049) 541-2121

### **EC Declaration of Conformity**

We,

Terumo (Philippines) Corporation  
124 East Main Avenue, Laguna Technopark  
Binan, Laguna, Philippines

whose single Authorized Representative:

Terumo Europe N.V  
Interleuvenlaan 40, 3001 Leuven, Belgium

Being the manufacturer, herewith declare that the products:

Terumo<sup>®</sup> Syringe with Needle  
Terumo<sup>®</sup> Syringe without Needle  
(with the attached product codes)

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class IIa according to Annex IX of the Directive 93/42/EEC. It bears the mark

**CE0197**

The product concerned has been designed and manufactured under a quality management system according to Annex V of Directive 93/42/EEC.

Compliance of the designated product with the Directive 93/42/EEC has been assessed and certified by the Notified Body

**TÜV Rheinland LGA Products GmbH**  
**Tillystraße 2, 90431, Nürnberg, Germany**

Certificate No.: DD 60108472 0001

Issue date: 2016 – 02 – 12

Expiry date: 2021 – 02 – 11

following the procedure relating to the "EC Declaration of Conformity" set out in Annex VII, combined with the provisions set out in Annex V "Production Quality Assurance" of Directive 93/42/EEC.



# Dossier d'information Euro Pharmat

## DISPOSITIF MEDICAL

This Declaration of conformity is valid in connection with the release document for the respective batch of produced devices.

This Declaration of Conformity covers all medical devices as specified in the product list belonging to this declaration and is only valid in connection with a batch specific Certificate of Compliance for all products concerned bearing the CE mark.

The above mentioned declaration of conformity is exclusively under the responsibility of

Terumo (Philippines) Corporation

Philippines, 02/17/14

**Place, date of issuance**

  
Alvin Robles

**Management Representative**

**Declaration of Conformity**  
**Terumo Syringe**  
**List of Product Codes**

**Terumo Syringe with Needle**

Volume mL	Product Description	Needle Size Gauge x Length	Product Codes	Lot Number
1	Syringe with Needle	25G x 5/8"	SS+01T2516	140116D
			SS+01T2516M	140213D
			SS+01T25161	131222D
			SS+01T25166	140203D
		26G x 3/8"	SS+01T2609	130907D
		26G x 1/2"	SS+01T2613	140207S
			SS+01T26131	130905D
		SS+01T26136	140208D	
	27G x 1/2"	SS+01T2713	140208D	
		SS+01T2713M	130713D	
Insulin Syringe	26G x 1/2"	SS+01H26131	130904D	
	25G x 5/8"	SS+01H25161	130802D	
2.5	Syringe with Needle – Luer Tip	21G x 5/8"	SS+02S21161	140918Y
		21G x 1"	SS+02S21251	140820Y
		21G x 1 1/2"	SS+02S21381	140916Y
		22G x 1 1/2"	SS+02S22381	140915Y
		23G x 1"	SS+02S23251	140708Y
3	Syringe with Needle – Lock Tip	20G x 1"	SS+03L2025M	100205F
		20G x 1 1/4"	SS+03L2032M	130511F
		20G x 1 1/2"	SS+03L2038M	130612F
		21G x 1"	SS+03L2125	121017F
			SS+03L2125M	131203F
		21G x 1 1/4"	SS+03L2132M	131226F
		21G x 1 1/2"	SS+03L2138	131203F
			SS+03L2138M	131203F
		SS+03L21386	141117F	
		22G x 1"	SS+03L2225M	130710F
		22G x 1 1/4"	SS+03L2232M	140209M
		22G x 1 1/2"	SS+03L2238	131206F
			SS+03L2238M	131217F
		23G x 1"	SS+03L2325	140212K
			SS+03L2325M	140213P
23G x 1 1/4"	SS+03L2332	140207P		
23G x 1 1/2"	SS+03L2338	121016F		
24G x 1"	SS+03L2425	131207F		
25G x 5/8"	SS+03L2516	140206F		



**Declaration of Conformity**  
**Terumo Syringe**  
**List of Product Codes**

Volume mL	Product Description	Needle Size Gauge x Length	Product Codes	Lot Number
			SS+03L2516M	130803F
3	Syringe with Needle – Luer Tip	21G x 1 ½"	SS+03S2138	131010A
			SS+03S2138M	131010A
		22G x ¾"	SS+03S2219	070831F
		22G x 1 ½"	SS+03S2238	130918A
		23G x 1"	SS+03S2325	140206A
		23G x 1 ¼"	SS+03S2332	131112F
		24G x 1"	SS+03S2425	130825A
		25G x 5/8"	SS+03S2516	131116A
5	Syringe with Needle – Lock Tip	20G x 1"	SS+05L2025M	091223C
		20G x 1 ¼"	SS+05L2032M	130801C
		21G x 1"	SS+05L2125	130923C
			SS+05L2125M	131018C
		21G x 1 ¼"	SS+05L2132M	140211R
		21G x 1 ½"	SS+05L2138	131016C
			SS+05L2138M	130615C
			SS+05L21386	141229C
		22G x 1"	SS+05L2225M	131120C
		22G x 1 ¼"	SS+05L2232	131130C
	SS+05L2232M		140210R	
	22G x 1 ½"	SS+05L2238	140207C	
		SS+05L2238M	130911C	
	23G x 1"	SS+05L2325	140130C	
	23G x 1 ¼"	SS+05L2332	140206C	
	Syringe with Needle – Luer Tip	21G x 1 ½"	SS+05S2138	130726C
			SS+05S2138M	131216C
SS+05S21381			140214C	
22G x 1 ¼"		SS+05S2232	110529C	
22G x 1 ½"		SS+05S2238	130810C	
		SS+05S22381	140216C	
23G x 1"		SS+05S2325	130811C	
23G x 1 ¼"	SS+05S2332	131113C		
	SS+05S23321	140221C		
10	Syringe with Needle – Lock Tip	20G x 1"	SS+10L2025	140213L
		20G x 1 ¼"	SS+10L2032M	130922L

Declaration of Conformity  
Terumo Syringe  
List of Product Codes

Volume mL	Product Description	Needle Size Gauge x Length	Product Codes	Lot Number
10		20G x 1 1/2"	SS+10L2038	061206E
		21G x 1"	SS+10L2125	131108E
	Syringe with Needle – Lock Tip	21G x 1 1/4"	SS+10L2132M	140214L
		21G x 1 1/2"	SS+10L2138	140123N
			SS+10L2138M	131204L
		22G x 1 1/4"	SS+10L2232M	140213L
		22G x 1 1/2"	SS+10L2238	131203L
			SS+10L2238M	131204L
	23G x 1"	SS+10L2325	140211L	
	Syringe with Needle – Luer Tip	20G x 1 1/2"	SS+10S20381	130607E
		21G x 1"	SS+10S2125	080104E
			SS+10S2138	140210E
		21G x 1 1/2"	SS+10S2138M	130702E
			SS+10S21381	130720E
		22G x 1 1/4"	SS+10S2232	121013E
		22G x 1 1/2"	SS+10S2238	131215E
SS+10S22381			130805E	
23G x 1"	SS+10S2325	140119E		

**Declaration of Conformity**  
**Terumo Syringe**  
**List of Product Codes**

**Terumo Syringe without Needle**

Volume mL	Product Description	Product Codes	Lot Number
1	Syringe without Needle	SS+01T	140204S
		SS+01TM	130206S
		SS+01T6	130731D
		SS+01T1	130815D
	Insulin Syringe without Needle	SS+01H1	130731D
		SS+01NA	131004S
2.5	Syringe without Needle – Luer Tip	SS+02S1	140215Y
3	Syringe without Needle – Lock Tip	SS+03L	130704K
		SS+03L1	140122K
		SS+03LM	131018A
		SS+03L6	140206P
	Syringe without Needle – Luer Tip	SS+03S	131012F
		SS+03S6	140206F
5	Syringe without Needle – Lock Tip	SS+05L	140127V
		SS+05L1	131203E
		SS+05L6	130918C
		SS+05LM	130731C
	Syringe without Needle – Luer Tip	SS+05S	130721V
		SS+05S1	140207V
10	Syringe without Needle – Lock Tip	SS+10L	140206N
		SS+10L1	131105N
		SS+10L6	131203N
	Syringe without Needle – Luer Tip	SS+10LM	140206N
		SS+10S	130621E
		SS+10S6	140207W
	Syringe without Needle – Eccentric Luer Tip	SS+10ES	140207E
		SS+10ESM	131009W
		SS+10ES1	140208E
20	Syringe without Needle – Lock Tip	SS+20L	131227B
		SS+20L1	140206B
		SS+20LM	140209B
	Syringe without Needle – Luer Tip	SS+20S	130513B
		SS+20ES	140206B
	Syringe without Needle – Eccentric Luer Tip	SS+20ES6	131130B
		SS+20ESM	140204B
		SS+20ES1	140203B

# Dossier d'information Euro Pharmat

## DISPOSITIF MEDICAL

Declaration of Conformity  
 Terumo Syringe  
 List of Product Codes

Volume mL	Product Description	Product Codes	Lot Number
30	Syringe without Needle – Lock Tip	SS+30L1	131105G
		SS-30L	150918G
50	Syringe without Needle – Eccentric Luer Tip	SS+50ES1	140206H
	Syringe without Needle – Catheter Tip	SS+50C1	140124H
	Syringe without Needle – Lock Tip	SS+50L1	140623H
60	Syringe without Needle – Catheter Tip	SS+60C	140107H
		SS+60CM	131026H
	Syringe without Needle – Lock Tip	SS+60L	150903H

No.DOC-PQB-TF-SS

Rev.08

### DECLARATION OF CONFORMITY

We, **TERUMO CORPORATION**  
44-1, 2-chome, Hatagaya, Shibuya-ku, Tokyo 151-0072, Japan

being the manufacturer of:

**TERUMO Syringe**

**Product : Hypodermic Syringe**

declare that the above products of **Class IIa** are in conformity with the provisions of the EC Council Directive 93/42/EEC of 14 June 1993, as amended, concerning medical devices, and have been subject to the conformity assessment procedure laid down in Article 11, 2 and 11, 3(a) of the Directive, relating to the "Full quality assurance" set out in Annex II, and by certification of Annex II, excluding Section 4 under the supervision of TÜV Rheinland LGA Products GmbH (Registration No.: HD 60077473 0001), Tillystraße 2, 90431 Nürnberg Germany, as Notified Body authorized by the German Competent Authority and carrying the Notified Body No. 0197.

Authorized European Representative :  
TERUMO EUROPE N.V.  
Interleuvenlaan 40, 3001 Leuven, Belgium

Object of the declaration: see appendix A

Tokyo, July 4, 2013  
(place and date of issue)

  
Hiroshi Nakagomi  
General Manager  
Quality Assurance Department  
TERUMO CORPORATION

 **TERUMO®**

No.DOC-PQB-TF-SS

Rev.08

### Appendix A - List of Code Number Structure

S S □ □ □ □ □ □ or  
1 2 3 4 5 7

S S □ □ □ □ □ or  
1 2 3 4 5

S S □ □ □ □ □ ○ ○ ○ ○ △  
1 2 3 4 5 6 7

1. Product series (product type) (one digit)  
SS : Syringe
2. Destination (Japan and overseas) (two digits)  
- : Japan \* : Overseas
3. Nominal capacity (product type) (two digits)  
02 : 2.5mL  
05 : 5mL  
10 : 10mL  
20 : 20mL  
30 : 30mL  
50 : 50mL
4. Cylinder head shape (one or two digits)  
S : Luer Slip tip  
L : Luer Lock tip  
ES : Eccentric Luer Slip tip  
C : Catheter tip
5. Others  
Z : Gamma sterilization  
E : Electron beam sterilization
6. Injection needle type (four digits)  
Upper two digits : Needle gauge  
Lower two digits : Needle length
7. Last digit  
1 : CE display

 TERUMO®

Page 2 of 2